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
DHR-0010 : TA Instruments Trios v3.1.4.3607

Experiment Instrument Engineering

Online Geometry Setup View Active Start Stop Go to Next Point Display Image View Log Controls Help

File Manager Instrument

Calibration

Inertia : 21.0048 $\mu\text{N.m.s}^2$

Last calibration date: 3/19/2014 2:15:43 PM

Calibration

Before calibration
Please ensure that no geometry is attached and that the spindle is free to rotate.

Calibration will take 30 seconds

Calibrate Accept Cancel

Rotational Mapping

Rotational mapping is on the geometry calibrations page Go To Geometry

Oscillatory Mapping

Current mappings

New mapping

Control panel

Idle

Name	Value	Units
Temperature	20.989	$^{\circ}\text{C}$
Set temperature		$^{\circ}\text{C}$
Torque	7.89255...	$\mu\text{N.m}$
Velocity	-0.0246...	rad/s
Displacement	386.161	rad
Axial force	0.312702	N
Viscosity	-1.9758...	Pa.s
Gap		mm

Gap

Closure profil: none

Set gap: mm

Apply

Gap

Environmental

Motor

Ready: 3:58:25 PM Gap monitor mode: none Smart swap : Smart swap read disabled Oven doors : partly open

calibration1 - ... DHR-0010 : TA I...

15:58 2014/4/28

TA

点击Calibrate

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DHR-0010 : TA Instruments Trios v3.1.4.3607

Experiment Instrument Engineering

Online Geometry Setup View Active Start Stop Go to Next Point Display Image View Log Controls Help

File Manager Instrument

Calibration

Instrument

Inertia : 21.0048 $\mu\text{N.m.s}^2$

Last calibration date: 3/19/2014 2:15:43 PM

Calibration

Calibration will take 30 seconds

Calibrating

Calibrate Accept Cancel

Rotational Mapping

Rotational mapping is on the geometry calibrations page Go To Geometry

Oscillatory Mapping

Current mappings

New mapping

Control panel

Calibrating

Name	Value	Units
Temperature	20.997	$^{\circ}\text{C}$
Set temperature		$^{\circ}\text{C}$
Torque	4.62316	$\mu\text{N.m}$
Velocity	-1.05874	rad/s
Displacement	387.088	rad
Axial force	0.312858	N
Viscosity	-2.6939...	Pa.s
Gap		mm

Gap

Closure profil: none

Set gap: mm

Apply

Gap

Environmental

Motor

Ready: 3:59:00 PM Gap monitor mode: none Smart swap : Smart swap read disabled Oven doors : partly open

calibration2 - ... DHR-0010 : TA I...

15:59 2014/4/28

TA

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DHR-0010 : TA Instruments Trios v3.1.4.3607

Experiment Instrument Engineering

Online Geometry Setup View Active Start Stop Go to Next Point Display Image View Log Controls Help

File Manager Instrument

Calibration

Instrument

Inertia : 21.0046 $\mu\text{N.m.s}^2$

Last calibration date: 4/28/2014 3:59:37 PM

Calibration

Calibration completed

New inertia 20.9978 $\mu\text{N.m.s}^2$

Calibrate Accept Cancel

Rotational Mapping

Rotational mapping is on the geometry calibrations page Go To Geometry

Oscillatory Mapping

Current mappings

New mapping

Control panel

Calibrating

Name	Value	Units
Temperature	21.005	$^{\circ}\text{C}$
Set temperature		$^{\circ}\text{C}$
Torque	0.256524	$\mu\text{N.m}$
Velocity	-0.800974	rad/s
Displacement	336.418	rad
Axial force	0.315249	N
Viscosity	-1.9758...	Pa.s
Gap		mm

Gap

Closure profil: none

Set gap: mm

Apply

Gap

Environmental


Motor

Ready: 4:00:11 PM Gap monitor mode: none Smart swap : Smart swap read disabled Oven doors : partly open

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

- 点击go to geometry或者点击calibrate—geometry calibrations 进入夹具校正窗口！

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DHR-0010 : TA Instruments Trios v3.1.4.3607

Experiment Instrument Engineering

Calibrate Rotational mapping Oscillatory mapping Options Load sample Service

File Manager

Geometries

- 25mm parallel plate, E...
- 25mm parallel plate, P...
- 40mm parallel plate, P...
- 40mm parallel plate, P...
- 60mm 1.999° cone plat...

25mm par...late Steel

25mm parallel plate, Peltier plate Steel

Inertia : 2.07589 $\mu\text{N.m.s}^2$ ✓

Last calibration date: 3/19/2014 2:16:42 PM

Calibration

Friction : 0.275919 $\mu\text{N.m/(rad/s)}$ ✓

Last calibration date: 3/19/2014 2:20:44 PM

Calibration

Gap Temperature Compensation : 0.7979 $\mu\text{m}/^\circ\text{C}$ ✓

Last calibration date: 3/19/2014 3:16:18 PM

Note: this calibration is only required for temperature ramps and temperature sweeps.

Calibration

Rotational Mapping ⚠

Last calibration date: 3/19/2014 2:33:36 PM

Calibration

Mapping settings

Bearing mapping type: standard

Number of iterations: standard, fast, precision

Mapping may take up to 5 minutes to complete.

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Environmental

Motor

Ready: 4:18:34 PM Gap monitor mode: none Smart swap : Smart swap read disabled Oven doors : fully open

16:18 2014/4/28

TA

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DHR-0010 : TA Instruments Trios v3.1.4.3807

Experiment Instrument Engineering

Calibrate Rotational mapping Oscillatory mapping Options Load sample Service

File Manager 25mm par...late Steel

Geometries

- 25mm parallel plate, E...
- 25mm parallel plate, P...
- 40mm parallel plate, P...
- 40mm parallel plate, P...
- 60mm 1.998° cone plat...

25mm parallel plate, Peltier plate Steel

Inertia : 2.07589 $\mu\text{N.m.s}^2$ ✓

Last calibration date: 3/19/2014 2:16:42 PM

Calibration

Friction : 0.275919 $\mu\text{N.m/(rad/s)}$ ✓

Last calibration date: 3/19/2014 2:20:44 PM

Calibration

Gap Temperature Compensation : 0.7979 $\mu\text{N}/^\circ\text{C}$ ✓

Last calibration date: 3/19/2014 3:16:18 PM

Note: this calibration is only required for temperature ramps and temperature sweeps.

Calibration

Rotational Mapping ⚠

Last calibration date: 3/19/2014 2:33:36 PM

Calibration

Mapping settings

Bearing mapping type: standard

Number of iterations: standard

Mapping may take up to 5 minutes to complete.

Options

☒ Display mapping prompt when the geometry changes

Calibrate Cancel

Control panel

Idle

Name	Value	Unit
Temperature	21.118	$^\circ\text{C}$
Set temperature		$^\circ\text{C}$
Torque	-0.719267	μN
Velocity	2.60681	rad
Displacement	660.107	rad
Axial force	-0.552268	N
Viscosity	-7.1948	Pa
Gap		μm

Gap

Closure profile: none

Set gap: Apply

Gap

Environmental

Motor


Ready: 4:18:34 PM Gap monitor mode: none Smart swap : Smart swap read disabled Oven doors : fully open

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


Discovery HR-1 : TA Instruments Trios

Home View Instrument Engineering

Paste HR-1 Experiment File Manager Display Document New Overlay Document New Analysis Document Edit New Geometry

File Manager [Experiment 1]

Experiments [Experiment 1]

Sample:

Name: Operator: Project: Notes:

Advanced

File Name: C:\Documents and Settings\All Users\Application Data\TA Instruments\TRIOS\Data\...

Geometry: 40mm parallel plate, Peltier plate Steel

Procedure:

Name:

1: Oscillation Frequency Control

Oscillation Flow Step (Transient) Other Conditioning

25 °C 0.0 s

Inherit set point Wait for temperature

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